

## Collapsed Z Stacks and 3D Projection with Ziess Images

For Images Acquired on the 510:

### Collapsed Z Stack:

1. Go to **Projection** in the main menu bar.
2. For Turning Axis select **Y**.
3. First Angle is set at **0**.
4. Number of projections is set at **1**.
5. Difference Angle does not matter what it is set at.
6. Click **Apply**.

Image will have the same name but with an \* added on the end of the file. You can change the name if you wish by clicking the save button in that window.

### Projection of Collapsed Z Stack Images:

1. Go to **Projection**.
2. In the **Transparency Tab** check for mode to be in **Maximum**.
3. In the Projection Tab select Y for the Turning Axis
4. For **Number of Projections and Difference Angle** the number must equal how much you wish to have the image turn around.

Example: If we want the image to turn 180 degrees we make the Number of projections = 18 and the Differentiation Angle = 10. For 360 degrees you make the Number of projections = 36 and the Differentiation angle = 10.

5. Check **Preview** if you wish to check how your image will turn out. You must use the slider to see the rotation.

6. Click **Apply**.

Again, image will have same name but with a \* added to the end of the file.

**If doing this on the microscopes you must first select 3D view in the menu to then select Projection.**

For Images Acquired on the 710:

Collapsed Z stack:

1. In the **Processing Tab** look for **Maximum Intensity Projection** [first one].
2. Click **Apply**.
3. Save image as a different name.